



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Operation Theatre Lights

Product Code :
SAMCOMEDICAL200001A001

Description :

Ceiling Operation Theatre Lights: Dome Diameter : 20in Light intensity : 180000 Lux Halogen Bulb : 24V - 150W
Reflector : Diachronic Colour Temp. : 4300K Temp. Increased at head area : 2~4C Finish : Oven Baked Ceiling
Operation Theatre Lights, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Samco Medical

Website: www.samcomedical.com,