



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Vaccine Carriers

**Product Code :**  
SAMCOMEDICAL220001A001

**Description :**

Vaccine Storage Capacity : 1.6 Litres Weight Fully Loaded : 4.2 Kgs. Weight Empty (with empty ice pack) : 2.2 Kgs. Insulation Thickness : 35/40 mm Approx. External Dimensions(WxDxH) : 25.0x25x30.0 cms Approx. Internal Dimensions(WxDxH) : 16.5x16.5x22.0 cms Vaccine Storage Dimension : 10.0x10.0x16.0 cms Number of Ice Packs supplied : 4 Volume of Ice Packs : 0.40 Litres Conforming to WHO Specifications : E4/VC.2 Cold Life Requirement as per WHO at 43 $\frac{1}{2}$ C  $\frac{1}{2}$  1 $\frac{1}{2}$ C without opening : Minimum 24 hrs. Actual Cold Life without opening at 43 $\frac{1}{2}$ C  $\frac{1}{2}$  1 $\frac{1}{2}$ C : 40 hrs.19 mins. Vaccine Storage Capacity : 0.90 Litres Vaccine Storage Capacity : 0.86 Litres Vaccine Carrier, Vaccine Carrier Exporters, Vaccine Carrier Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),