



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Suspended OT Lights

Product Code :
SAMCOMEDICAL4008A004



Description :

Ceiling Suspended OT Lights Ceiling Suspended OT Lights are 19 $\frac{1}{2}$ dia metallic dome fitted with halogen bulb (1x 24V/150W) & heat absorbing colour correcting filters Ceiling Suspended OT Lights are Provided with plastic diffusers for shadow reduction. Spring balanced Ceiling Suspended OT Lights are Provided with sterilizable handle and low voltage supply unit Ceiling Suspended OT Lights are Light Intensity: 25,000 Lux + 10% Ceiling Suspended OT Lights, LED OT Light exporters, LED OT Light manufacturers



Samco Medical

Website: www.samcomedical.com,