



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Binocular Stereo Microscope

**Product Code :**  
SAMCOMEDICAL4011A002



**Technical Specification :**

Binocular Stereo Microscope The Binocular Stereo Microscope is designed for younger students of elementary education. It meets the needs of Science teaching and demonstration to beginners at school level. Being light weight & handy, it helps carrying around and is simple to use, yet having features of a compound Microscope Precision Optical Glass Lenses and metal components are engineered for years of reliable performance. Observation Head - Binocular head inclined at 45-deg. with adjustable inter-pupillary distance. Focussing - By Rack and Pinion mechanism. Base - Fitted on transillumination base suitable for artificial and transmitted arrangement and elbow hand rest. Magnification - Standard magnification 20X & 40X or 10x & 30x. Objectives - 2 sets of objectives 2X and 4X, selected by sliding nosepiece. Eyepieces - Standard extra wide field eyepieces WF-10X or 15X are provided. Binocular Stereo Microscope AEI-A, Binocular Loupes exporters, Binocular Loupes manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),