



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dental Surgical Microscope

Product Code :
SAMCOMEDICAL4012A002



Description :

Dental Surgical Microscope Owing to our experience, we have been successful in catering to the requirements of clients by offering quality Dental Surgical Microscope. All the products are widely acclaimed by the customers for their excellent performance, easy installation, low maintenance and longer service life. Apart from this, we provide customization as per the specifications laid down by the customers. Specification: Binocular Tube: 45° 1/2, F=170mm Magnification: 0.4x, 0.63x, 1.0x, 1.6x & 2.5x with five step magnifier Eye Piece: WF 10x/16mm with eye guards Working Distance: F=200mm objective lens Focusing: Motorized fine focusing Light Source: 15v, 150w twin halogen reflector lamp Built in Filters: Green & Cobalt Blue Vertical Movement of Arm: 600mm Microscope Carriers: XYZ Coupling or Z-axis Ophthalmic carrier Optional: Assisto scope, Xpert Vu Camera Module, Double Beam Splitter, 0-210° 1/2 tiltable ergo tube, eye piece WF 12.5x/16mm Dental Surgical Microscope, Dental Microscopes exporters, Dental Microscopes manufacturers



Samco Medical

Website: www.samcomedical.com,