



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dental Surgery Microscopes

Product Code :
SAMCOMEDICAL4012A003



Description :

Dental Surgery Microscopes Innovative designed, highly versatile with latest 45° inclined optical system 15x wide field eye pieces 3 Step magnification (5x, 10x & 20x) or 5 Step magnification (0. 4x, 0. 63x, 1x, 1. 6x & 2. 5x) or Continuous zoom (5x ~ 25x) Beam splitter with 1/2° mount CCD digital camera Working distance : f= 200mm objective lens Interpupillary distance adjustment: 55mm to 75mm Focusing 1/2° manual or motorized with foot control LCD monitor/screen : mount on microscope to display live surgery Illumination section: Light source : 15v-150w / 24v-250w twin halogen reflector lamp Illumination control : continuous adjustment Co-axial through-the-lens fiber optic light guide, white bright light Field of illumination : 70mm Filter : cobalt blue, red-free, heat absorbing filter Stand section: Type: mobile floor stand on four caster wheel for easier handling and absolute stability Vertical height adjustment range of counter balanced pantographic arm: 300 to 560mm Rotation of arm : 320 with lock General information: Base size: 600mm x 550mm Weight: 65 kgs. Approx. Height: 1550mm Dental Surgery Microscopes, Dental Microscopes exporters, Dental Microscopes manufacturers



Samco Medical

Website: www.samcomedical.com,