



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Portable Surgical Microscopes

Product Code :
SAMCOMEDICAL4013A001

Description :

Portable Surgical Microscope Backed by our efficient workforce, we are able to design and develop a wide range of Portable Surgical Microscope. These are developed with 45 deg inclined optical system, 3step magnification optical head (5x, 10x, 20x), 15x wf eye piece, working distance f = 200 mm (300 mm, 400 mm) objective lens, 15v/150w or 24v/250w cold light source with continuous adjustments illumination low & high intensity control, manual fine focusing. Our products are mainly used for dental surgery, plastic surgery & neuro surgery.
Features: ½ Lightweight ½ Easy handling ½ Low maintenance Portable Surgical Microscope, Portable Surgical Microscope exporters, Portable Surgical Microscope manufacturers



Samco Medical

Website: www.samcomedical.com,