



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
General Surgery Microscopes

**Product Code :**  
SAMCOMEDICAL4014A001

#### Description :

General Surgery Microscope We provide our customers with General Surgery Microscope that is available in different models and sizes at most competitive prices. These products are designed in compliance with international quality norms & standards, using high-grade raw material and latest technology. Our products are mainly used in various surgeries like neuron- surgery, gynecology, plastic surgery and general surgery.

Features: Durability Reliability Dimensionally accurate Accurate testing results Technical specifications:

Microscope section: Innovative designed, highly versatile with latest 45° 1/2 inclined optical system 15x wide field eye pieces 3 Step magnification ( 5x, 10x & 20x) or 5 Step magnification (0. 4x, 0. 63x, 1x, 1. 6x & 2. 5x) or Continuous zoom (5x ~ 25x) Beam splitter with 1/2 inch mount CCD digital camera Working distance: f= 200mm objective lens Interpupillary distance adjustment: 55mm to 75mm Focusing 1/2 inch manual or motorized with foot control LCD monitor/screen: mount on microscope to display live surgery Illumination section: Light source : 15v-150w / 24v-250w twin halogen reflector lamp Illumination control : continuous adjustment Co-axial through-the-lens fiber optic light guide, white bright light Field of illumination : 70mm Filter : cobalt blue, red-free, heat absorbing filter Stand section: Type: mobile floor stand on four castor wheel for easier handling and absolute stability Vertical height adjustment range of counter balanced pantographic arm: 300 to 560mm Rotation of arm : 320 with lock General information: Base size : 600mm x 550mm Weight : 65 kgs. Approx. Height : 1550mm Electrical power : AC 220V, 240V, 50/60hz or 110 V on request Power consumption : 220va or 110 V Optional accessories: Monocular assisto scope Objective lens 250mm, 300mm & 400mm General Surgery Microscope, General Surgery Microscope exporters, General Surgery Microscope manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),