



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Colposcope 3 Step Magnification

Product Code :
SAMCOMEDICAL4015A001

Description :

Colposcope 3 Step Magnification Our Colposcope 3 Step Magnification is commonly used by Gynaecologists for directly visualization of vulva, vagina and cervix in addition to other areas. It is very flexible and can adjust the binocular eye pieces to get a fine focus. Features Objective Lens Moveable Binocular Eye Pieces Height is Adjustable Fine Focusing Tilting Arrangement HighLights High Precision Optics Reliable for In-office Surgical Application. Also Excellent for Operation Theatre and Lab Work. Very Economically priced. Specifications Colposcope Tilting Adjustment: 30o Upward, 30o Downward Working Range (Fine Focusing): 50mm Lamp: 15V / 150W, 24V / 250W Total All-Inclusive Weight: 22 Kg. Filter: Red Free Green Individual Components Eye Piece: 15x, 12.5x Objective Lens: 200mm, 300mm, 400mm Field Of Vision: 32mm Magnification: 5x, 10x, 20x - 3 step Ross Slide: 50mm Light Source: Halogen Fibre Optic.Single or Double Outlet Optional Accessories Beam Splitter With CCD Camera Adapter Table Clamp Attachment CCD Camera Optional accessories are listed as separate items in "Accessories" section under products, there can be bought separately. Applications Gynaecological (Including Vaginal and Uterine) Diagnosis and/or Surgery or Post-op. Inspection. Detection of Uterine Cancer, Vaginal or Uterine Erosions, Tumors, Congestion of Inflammation. Very Accurate Biopsy of Uterus or Vagina, (With Help Of Clearly Magnified View.)



Samco Medical

Website: www.samcomedical.com,