



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Colposcope 5 Step Magnification

**Product Code :**  
SAMCOMEDICAL4015A003



**Description :**

Colposcope 5 Step Magnification Colposcope 5 Step Magnification provides a smooth zoom magnification continuous fashion. It provides a great degree of focus and resolution. It is very flexibility and can adjust the binocular eye pieces to get a fine focus. Features High Lights Objective Lens High Precision Optics Moveable Binocular Eye Pieces Reliable for In-office Surgical Application. Height is Adjustable Also Excellent for Operation Theater and Lab Work. Fine Focusing Very Economically priced. Tilting Arrangement Specifications Colposcope Tilting Adjustment : 30°½ Upward , 30°½ Downward Working Range (Fine Focusing) : 50mm Lamp : 15V / 150W, 24V / 250W Total All-Inclusive Weight : 22 Kg. Filter : Red Free Green Individual Components Eye Piece : 15x, 12.5x Objective Lens : 200mm, 300mm, 400mm Field Of Vision : 32mm Magnification : 0.5x, 7.5x, 12.5x, 20x, 30x- 5 step Cross Slide : 50mm Light Source : Halogen Fiber Optic.Single or Double Outlet Optional Accessories Beam Splitter With CCD Camera Adapter Table Clamp Attachment CCD Camera Optional accessories are listed as separate items in "Accessories" section under products, there can be bought separately. Applications Gynecological (Including Vaginal and Uterine) Diagnosis and/or Surgery or Post-op. Inspection. Detection of Uterine Cancer, Vaginal or Uterine Erosions, Tumors, Congestion of Inflammation. Very Accurate Biopsy of Uterus or Vagina, (With Help of Clearly Magnified View.) Colposcope 5 Step Magnification, Gynecology Equipment exporters, Gynecology Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

