



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fiber Optic Otoscope

Product Code :
SAMCOMEDICAL4015A004



Description :

Fiber Optic Otoscope Otoscope with Fiber Optic (F.O.) illumination constructed of light, high-impact Polycarbonate and 3x magnification scratch resistant acrylic viewing window. Output for pneumatic membrane test. Otoscope Set complete with 10 mini Speculum Size 2.5mm and 4mm. Set complete with C size Dry battery Handle. Fiber Optic Otoscope, Gynecology Equipment exporters, Gynecology Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,