



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Digital Colposcope

Product Code :
SAMCOMEDICAL4016A001

Description :

Digital Colposcope We are engaged in offering our clients with Digital Colposcope that is available in different models and sizes at market leading prices. These are designed in compliance with international quality norms & standards, using high-grade raw material and latest technology. Apart from this, we check all the products on certain parameters to ensure flawlessness. Digital Colposcope: High-person optics Reliance for in-office surgical application Also excellent for operation theater and lab work Very economically priced Also conveniently moveable for one place to another, even by lightly-built/lady gynaecologists Applications: Gynecological (Including vaginal and uterine) diagnosis Detection of uterine cancer,vaginal or uterine erosion, tumors, congestion of inflammation. Very accurate Biopsy of uterus or vagina Specifications: Tilting Adjustment 30 Upward, 30 Downward Working Range(Fine Focusing 50mm Lamp 15V/150W,24V/250W Total all-inclusive weight 22kg. filter Red Free Green Individual Components: Eye Piece 15 x,12.5x Objective Lens 200mm,300mm,400mm Field of vision 32mm Magnification 5x,10x,20x-3 step Model maxcp 1 5x7.5x,12.5x,30x-5 Step Model-VICS-30 5x to 25x continuous Zoom Model-maxcp-1 Cross Slide 50mm Light Source Halogen Fibre, Optic, Single Or Double outlet Optional Accessories Beam Splitter With CCD Camera Attachment Digital Colposcope, Gynae Colposcope exporters, Gynae Colposcope manufacturers



Samco Medical

Website: www.samcomedical.com,