



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Slit Lamp Zoom Mag

Product Code :
SAMCOMEDICAL4002A004



Description :

Slit Lamp Zoom Mag Magnification: 5x,10x,15x,20x & 25x continuous zoom Inter pupillary Distance: 52mm to 80mm Camera: CCD Digital Camera Illumination Section: Slit Width Adjustment: 0 to 12mm without steps Slit Length Adjustment: 0 to 12mm continuous Diameter of Diaphragm: 12mm, 10mm, 8mm, 5mm, 3mm, 2mm, 1mm and 0.2mm Angle of Rotation of Slit: $1\frac{1}{2}$ to $90\frac{1}{2}$ Tilt of Slit: Horizontal in steps of $0\frac{1}{2}$, $5\frac{1}{2}$, $10\frac{1}{2}$, $15\frac{1}{2}$, Vertical in steps $0\frac{1}{2}$, $5\frac{1}{2}$, $10\frac{1}{2}$, $15\frac{1}{2}$, $20\frac{1}{2}$ Filters: Cobalt Blue, Green, Neutral Density, Heat absorbing Light Source: Halogen Lamp 12v-50w Intensity Control of Illumination: Low, Medium and High



Samco Medical

Website: www.samcomedical.com,