



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Slit Lamp Z Excel

**Product Code :**  
SAMCOMEDICAL4002A005



**Description :**

Slit Lamp Z Excel Slit Lamp Excel is often used as part of a routine eye exam. It is used to detect disorders in the front of the eye, such as clouding of the lens (cataracts), glaucoma, or abnormalities of the cornea. It provides a magnified, three-dimensional view of the structures within the eye. Specifications:- Binocular Microscope Eye Pieces : 10x, 15x Objective : 1x, 1.6x Total Magnification : 10x, 16x and 24x Illumination Unit Slit Image Rotation : 0 to 180 Filter Disc : Cobalt Filter, Green Filter, Yellow Filter Natural Density and Open Aperture Slit Disphragm Disc : Wedge Shaped Diaphragm Of Infinitely Variable Slit Lengths Halogen Lamp : 6 Volt 20 Watts Standard Accessories 15x Eyepieces Plastic Breath Shield Dust Cover and Fuse Testing Rod Replacing Lamp One Piece For Illumination Unit Optional Accessories:- Beam Splitter With CCD Camera Adaptor CCD Camera Power Table / Mechanical Table Optional accessories are listed as separate items in "Accessories" section under products, there brought separately. Slit Lamp Z Excel, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),