



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Slit Lamp Max

Product Code :
SAMCOMEDICAL4002A007



Description :

Slit Lamp Max Slit Lamp Max features are computer designed, multi-coated optics that provide unsurpassed viewing of ocular structures, the operator's field of view is maximized for precise visualization of ocular defects at the periphery. Specifications Binocular Microscope Illumination Unit Slit Image Rotation : 0 to 180° ½ Tilting Illumination : Up to 20° ½ Natural Density and Open Aperture Slit Disphragm Disc : Wedge Shaped Diaphragm Of Infinitely Variable Slit Lengths Halogen Lamp : 6 Volt 20 Watts Standard Accessories 15x Eyepieces Plastic Breath Shield Replacing Mirror Hruby Lens Dust Cover and Fuse Testing Rod Replacing Lamp one piece for illumination unit Optional Accessories CCD Camera Beam Splitter With CCD Camera Adaptor Power Table / Mechanical Table Optional accessories are listed as separate items in "Accessories" section under products, there brought separately. Specification :- Slit Lamp Max Eye Pieces 10X, 15X Objective 1X, 1.6X Magnification 10X, 16X & 24X Filter Disc Cobalt Filter, Green Filter Yellow Filter Slit Lamp Max, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,