



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Diagnostic Sets

**Product Code :**  
SAMCOMEDICAL4002A013



**Description :**

Diagnostic Sets Diagnostic Set Includes: - One Ophthalmoscope Head with 6 different apertures micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter. 28 focusing lenses with a range of -25 to +40 diopters. Rubber brow rest prevents scratching of eyeglasses. Otoscope Head with 3x magnification lens fitted with 3.5v Halogen Lamp. One Halogen Otoscope Head with Reusable Ear Specula 3mm & 4mm, One Halogen Ophthalmoscope Head with 6 Apertures, One Spreadable Nasal Speculum, One Bent Arm Illuminator, 2 Laryngeal Mirrors 4 No. & 5 No., One Tongue Depressor, One Tongue Depressor Holder, One "C" cell Dry Battery Handle. Set complete in a Case. Diagnostic Sets, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),