



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ophthalmoscope Standard

Product Code :
SAMCOMEDICAL4002A024



Description :

Ophthalmoscope Standard We bring in advanced quality standard ophthalmoscope for sale. It is widely used in various diagnostic and surgical applications of the eyes. The standard ophthalmoscope head we make available, uses high quality halogen light for long lasting illumination. Our medical ophthalmoscope is remarkable for viewing the true color of eyes tissues. They are accessible in six different apertures i.e. micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter. In addition to it, our standard ophthalmoscopes come with rubber brow that prevents any scratch on the eyeglasses. Further, we are among the leading suppliers of standard ophthalmoscope. The Standard Ophthalmoscope uses halogen light for true tissue color and consistent, long-lasting illumination. It is versatile and reliable for general use. Features Halogen light for true tissue color and consistent, long-lasting illumination. 6 apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter. 24 focusing lenses with a range of -25 to +40 diopters. Rubber brow rest prevents scratching of eyeglasses. Ophthalmoscope Standard, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,