



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Otoscope

**Product Code :**  
SAMCOMEDICAL4002A025



**Description :**

Otoscope 2.5v Halogen Otoscope with Metal Head in Mat finishes with Pin contact fitting. Optical glass viewing window with 3x magnification. Movable Swiveling viewing window, permanently-attached. Metal reusable speculum size 2, 3, 4mm. Output port for pneumatic test. Rheostat to adjust the intensity of light. C size Battery Handle. Available in Wooden Case or in zipper pouch Otoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),