



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Streak Retinoscope

**Product Code :**  
SAMCOMEDICAL4002A029



**Description :**

Streak Retinoscope is exceptionally bright. It is a Streak and Spot Retinoscope in one and is converted by simply changing the bulb. Retinoscope with D size battery handles in case or zipper pouch packing. Rotation and Vergences movements independently operated for those who prefer separate controls. Continuous 360° 1/2 rotation for easy and accurate positioning of streak for axis determination. Streak Retinoscope Streak Retinoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),