



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ophthalmometer Instrument

**Product Code :**  
SAMCOMEDICAL4002A032



**Description :**

Ophthalmometer Instrument Our business group is one of the prominent Ophthalmometer Instrument manufacturer and supplier. Manufactured using quality basic material, these ensures longer service life and reliable performance. Astigmatism is a condition of the eye in which the cornea or lens is misshapen and can result in vision problems. Specification: Measuring range Dioptres from 30,00 to 60,00 D - radius from 5,6 to 11,3 mm Resolution Dioptres 0,25 D - radius 0,05 mm Useful measuring diameter 3,4 mm equal to 9 square mm Rating scale Table & International Lines of sight According to Javal Instrument base Cross-wise movement to be stopped Chin rest Adjustable height with shutter device Feeding voltage 120 Vac /60 Hz - 230 Vac/50 Hz Instrument voltage 6 VaC Bulbs of the lines of sight 6 V - 10 W - SV 8,5 - 8 Bulb of the illuminating rating scale 6 V - 3 W - E 10 Fuses 5x20 mm 120 V - 250 mA T230 V - 125 mA T Absorbing power (VA) 24 VA Ophthalmometer Instrument, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),