

Product Name :
Refraction Unit

Product Code :
SAMCOMEDICAL4002A037



Description :

Refraction Unit This Refraction Unit has a compact design which ensures optimal functionality of the instruments and can be installed in a minimum space. The refraction unit is so planned that it simplifies all functions offering the best cost performance ratio having all quality features of a plan tech products. Features: ½ Fully Automatic Operation ½ CE Certified Noiseless German Motors for up, down & Reclining ½ Easy to Operate ½ Excellent quality & well Rounded design ½ Long life ½ Adjustable foot rest ½ Optimal Utilization of space Specifications: Elevation Strokes 200 mm Seat Height 480mm to 580mm Inclination of Back Rest 90 to 180 Degree Rotation (Optional) 180 Degree Power Source 220 V A/c 50Hz Power Consumption 660 VA Power Supply 0, 1. 5, 2. 5, 3, 4, 6, 7. 5, 9 & 12V Sliding Tray Loaded Weight 30 Kg on right & left side Height from the floor 900mm Table Top Size 810mm x 405mm Refraction Unit, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,