



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Riester Schiotz C Tonometer

Product Code :
SAMCOMEDICAL4002A041



Description :

Riester Schiotz C Tonometer High quality agate bearing for an extremely long service life Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions Perfect reading of the scale with red pointer All vital parts are made of stainless steel, other parts are chrome-plated Complete in a black deluxe case with velvet-look inserts Riester Schiotz C Tonometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,