



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Portable Slit Lamp

Product Code :
SAMCOMEDICAL4002A043



Description :

Portable Slit Lamp) latest innovation 2) skillful manufacture 3) quality and quantity assured 4) Handful operation 5) Durable using Eye piece 10X,16X Total Magnification 10X,16X Adjustable Diopter $\pm 7D$ Work distance 80mm Slit Width(mm) 0~10mm,Continuously adjustable (at 10mm, slit becomes a circle) Slit Angle Horizontal circumference $\pm 30^\circ$ Filter Green, cobalt blue, color temperature transformation angle Illumination bulb White LED Voltage & Frequency: AC 100~240V 60~50Hz Power 18VA Portable Slit Lamp, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,