



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Semi Automatic Capsule Filling Machine

**Product Code :**  
SAMCOMEDICAL4003A002



**Description :**

Semi Automatic Capsule Filling Machine We hold expertise in producing a range of semi automatic capsule filling machines that are made using superior quality raw material. Designed with precision, our range of machines are known to provide optimum automation level while filling capsules of all sizes with powder, granules and pellets. Further, we also make sure that these machines are quality tested to meet the international industrial standards. Salient features : PLC based machine A special design for change over of different size capsules with in 20 minutes. Fill weight accuracy +/-3%. Simple operation, operating training can be done in 15 minutes Low noise & maintenance. CGMP with documentation. Capsule Size:00,0,1&2,3 to Semi Automatic Capsule Filling Machine, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),