



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Tube Filling Machine

Product Code :
SAMCOMEDICAL4003A003



Description :

Tube Filling Machine We are offering tube filling machine which is used to seal and fill a wide range of soft plastic tubes, metal tubes and flexible compound tubes. The machines manufactured by us are well known to provide reliable, flexible and user friendly operations. The machine is designed in such a way that provides good access for maintenance and give operator excellent visibility at all the major machine functions. Compactly constructed, these machines are equipped with high accuracy filling systems that provides clean and accurate filling at optimum speed providing high quality products. Features : Automatic Tube Loading System Auto Tube Orientation No Tube No Fill System Folding & Coading at Individual Station Simple Mechanism & Operator Friendly Machine Combo Model Suitable For Both Aluminum and Lami/Plastic Tubes Out Put 70 to 80 Tubes Min - For Double Head M/c Out Put 35 to 40 Tubes Min For Single Head M/c Technical Specification : Filling Range 3 Gm To 250 Gm With Change Over. Tube Dia Upto 25 mm - For Double Head. Tube Dia Upto 45 mm - For Single Head. Tube Filling Machine, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,