



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Octagonal Blender

**Product Code :**  
SAMCOMEDICAL4003A005



**Description :**

Octagonal Blender Octagonal Blender is widely used in pharmaceutical industry, which demands gentle blending of dry granules-powder. Its structure is octagon shaped with a rectangular central portion and top & bottom polygonal frustum. Its salient features are as follows: Polygon shaped body. Contact parts in stainless steel 304 or 316. Manhole with gasket to fix baffles. Charging and discharging is through a butterfly valve.. Easy for cleaning in places. Lower power consumption. All moving parts are enclosed in guards. A safety railing with a limit switch for operator safety. Octagonal Blender, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),