



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Zero Hold Up Filter Press

Product Code :
SAMCOMEDICAL4003A006



Description :

Zero Hold Up Filter Press We are manufacturers & offering zero hold up filter press which is used for crystal clear filtration in pharmaceuticals, chemicals and other industries. Manufactured using the high grade components, these presses are dimensionally accurate and requires low maintenance. Designed in adherence with international standards, these machines are easy to operate and install. Available in different size and specifications, these are also customized as per specific requirements of the clients. Salient features : Sizes available 8, 14, 18, 24, & 33, diameter Moc ss 316 / 316L. Available models 8, 14, 18, 24, & 33, plate diameter. Available with centrifugal / gear / lobe pump. Fitted with pressure gauge, safety valve. Zero Hold Up Filter Press, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,