



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Rapid Mixer Granulator

**Product Code :**  
SAMCOMEDICAL4003A007



**Description :**

Rapid Mixer Granulator We present a versatile range of rapid mixer granulators which are used to mix and granulize the various ingredients in pharma and other allied industries. These machines provide better and fast mixing which leads to fast tableting speeds with least rejections and improved quality. The granulators are known for their high productivity with low operating cost. In addition, we also provides customized specifications as required by our client. Salient features : Models available from 3 L to 800 L working capacity Fully automatic operation based on PLC control. Better productivity & lesser mixing time. Contact parts are SS 316 / 316L Uniform mixing & granulation Less handling. More hygienic. CGMP with documentation. Rapid Mixer Granulator, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),