



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
POWDER AND MASS MIXER

Product Code :
SAMCOMEDICAL4003A008



Description :

POWDER AND MASS MIXER Suitable to mix dry and wet material. The mixer has all contact parts made out of 304/316 quality stainless steel with silicon rings and asbestos rope at both the ends, to prevent outside lubricants from entering into mixing chamber. The paddles on the shaft are so arranged as to give thorough and uniform mixing. The unit is supplied with tilting arrangement, reversible switch and safety lock-switch to stop the machine as soon as the cover is opened. The drive is through a suitable motor and gear box covered with S.S. The unit is as per GMP Standard. All the M.S. part are cladded with S.S. POWDER AND MASS MIXER, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,