



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
MULTI MILL

Product Code :
SAMCOMEDICAL4003A009



Description :

MULTI MILL The Multi Mill is useful for various purposes like granulation, shredding, pulverization etc. The unit consists of vertically mounted 3 HP motor which drives the motor shaft through bearing housing. The rotor has 12 Nos. swinging type knife blades and 2 Nos. scrapper blades. The- screen of suitable mesh is mounted outside the rotor and material is fed from S.S. Hopper. All the contact parts are made out of S.S. 304/316q. The unit is mounted on suitable M.S. Stand with castor-wheels and is provided with DOL starter, reversible switch and five numbers of Stainless Steel screens. (Size of Screens to be specified while placing the order). Same as above but as per GMP Std. M.S. Stand and all visible M.S. Parts are clad with S.S. Extra price for spare screen of any one size of (i.e. 12 / 10" 8/ 6/ 5/ 4/ 3/ 2/ 0.5mm) MULTI MILL, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,