



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
STAINLESS STEEL COATING PANS

Product Code :
SAMCOMEDICAL4003A011



Description :

STAINLESS STEEL COATING PANS Coating Pans are made out of stainless steel and are fitted on a drive arrangement having suitable oil filled worm gear box and motor. The unit is also equipped with suitable hot air blower driven by means of suitable motor with digital Thermostatic Control. The whole unit will be cladded with S.S. including flexible pipe of Blower is also made out of stainless steel. No painted parts visible as per GMP Std. STAINLESS STEEL COATING PANS, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,