



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
VACUUM TRAY DRYER

Product Code :
SAMCOMEDICAL4003A013



Description :

VACUUM TRAY DRYER

- Overall body of thick S.S. Plate with optional SS 304Q or SS 316Q.
- Option of hot water circulation by limped coil for outer body to avoid dew/moisture formation on the internal chamber of dryer.
- Option of providing receiver condenser or vacuum pump to recover the harmful chemicals or expensive solvents
- Maximum contact surface for efficient heat transfer to product by giving flattened bottom of the trays & top plate of heating shelves.
- Longer & durable heating pad plates made from thick SS sheet fully argon arc gas Welded & provided with S.S. 304Q square rod stiffeners
- Hinged 'C' Clamp with knobs all across the door to avoid leakages

VACUUM TRAY DRYER, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,