



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
DRUM MIXERS CAPACITY 200 LTRS.

Product Code :
SAMCOMEDICAL4003A014



Description :

DRUM MIXERS CAPACITY 200 LTRS Complete with stainless steel drum of 620mm dia x 650mm height made out of 304 quality 16 gauge stainless steel. The drum is complete with air tight lid with quick closing and opening clamps. The unit is driven by 1 HP, 1440 RPM Motor through 1 3/4" gear box DRUM MIXERS CAPACITY 200 LTRS, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,