



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
STRIP DEFOILING MACHINE

Product Code :
SAMCOMEDICAL4003A019



Description :

STRIP DEFOILING MACHINE This machine is suitable for defoiling the tablets/ capsules from the rejected foils to recover tablets, capsules as well as aluminium foil. The foils are fed in the machine and the tablets/ capsules are separated. The machine is provided with adjustable slitters suitable for strip width of 50 to 90mm and is supplied complete with stainless steel tablet collecting tray etc. All the shaft and rollers are made out of stainless steel and the unit is driven by a HP (3 Phase A/C 440 Volts) Motor STRIP DEFOILING MACHINE, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,