



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
STAINLESS STBEL CYLINDRICAL TANKS

Product Code :
SAMCOMEDICAL4003A021



Description :

STAINLESS STBEL CYLINDRICAL TANKS Made out of stainless steel (304/316 q.) 16 gauge with all the welding done by argon arc process. The tanks are polished and buffed from inside and outside to give mirror finish and are supplied complete with lid, one outlet valve with hose nipple, S.S. claddcd M.S. Stand having castor wheels etc. Model T-50 and T-100 are supplied without valves and stand. The dimensions are approximate and are subject to change depending on availability of S.S. Sheets STAINLESS STBEL CYLINDRICAL TANKS, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,