



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
RECACTION VESSEL

**Product Code :**  
SAMCOMEDICAL4003A022



**Description :**

RECACTION VESSEL Stainless Steel 304/316L quality.  $\frac{1}{2}$  MS/SS Jacket or limpet coil for steam heating / cooling.  $\frac{1}{2}$  Reaction Vessel can be provided with different types of Agitator like Anchor type / Propeller type / Paddle type Pitched blade turbine type of with high speed homogenizer.  $\frac{1}{2}$  Shaft entry sealed by Oil seal / Plain or Water cooled Stuffing Box / Single or Double Mechanical seal.  $\frac{1}{2}$  High density glass wool insulation & S.S. Cladding can be done on demand.. RECACTION VESSEL, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),