



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi Automatic Cap Sealing Machine Versatile Ropp Cap Sealing Machine

Product Code :
SAMCOMEDICAL4003A028



Description :

Semi Automatic Cap Sealing Machine Versatile Ropp Cap Sealing Machine

Technical Specification :

SEMI AUTOMATIC CAP SEALING MACHINE VERSATILE R.O.P.P. CAP SEALING MACHINE The unit is covered by S.S, sheets with motor and clutch arrangement to give positive sealing. Output 30 to 40 bottles/ minute. SEMI AUTOMATIC CAP SEALING MACHINE VERSATILE R.O.P.P. CAP SEALING MACHINE, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,