



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Surgical Microscope Direct Halogen

Product Code :
SAMCOMEDICAL4004A002



Description :

Surgical Microscope Direct Halogen Surgical Microscope Direct Halogen is also same as Surgical Microscope Cold Light but with the difference that the additional illuminating beam Halogen Lamp and a measurement beam which makes the patients surgery very clear. Specifications Microscope Section Binocular Tube : 90° 1/2 ENT & General Surgery Eyepieces : 12.5x/ 10x/ 16mm Dioptic Adjustment Magnification Changes : Built in Manual 3-Step Changer Total Magnification : 7.5x, 12.5x, 20x Field of View : 59mm, 37.8mm, 23.5mm, 14.5mm, 9.5mm Objective Lens : F = 200, 250, 300, 400mm (For General Micro-Surgical Fields) Working Distance : 200mm Interpupillary Distance Adjustment : 55mm to 75mm Illumination Section Type : Direct Halogen Light Source : 12V/50Watts Halogen lamp Illumination Control : Continuous Adjustment Field of illumination : 70mm Filter : Cobalt Blue, Red Free, Heat Absorbing Stand Section Type : Table Clamp Vertical Movement Of Counterbalanced Arm : 400 mm General Information Weight : 50Kg. Height : 1550mm Electrical Power : AC 220V, 240V, 50/60 Hz Power Consumption : 220 VA/110V Optional Accessories Beam Splitter With CCD Camera Attachment CCD Camera Optional accessories are listed as separate items in "Accessories" section under products, there brought separately. Surgical Microscope Direct Halogen, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,