



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ent microscope

Product Code :
SAMCOMEDICAL4004A003



Description :

Ent microscope Surgical Microscope Cold Light is a highly versatile with a Zooming Optical System which helps in focusing. This microscope has a high magnification for making the view of a particular part clear for surgery. Specifications Microscope Section Binocular Tube : 90° 1/2 ENT & General Surgery Eyepieces : 12.5x/ 10x/ 16mm Dipprtric Adjustment Magnification Changes : Built in Manual 5-Step Changer Total Magnification : 4x, 6x, 10x, 16x, 25x, 1:6 Field of View : 59mm, 37.8mm, 23.5mm, 14.5mm, 9.5mm Objective Lens : F = 200, 250, 300, 400mm (For General Micro-Surgical Fields) Working Distance : 200mm Interpupillary Distance Adjustment : 55mm to 75mm Illumination Section Type : By Optical Light Guide System Light Source : 15V/150Watts Halogen lamp Illumination Control : Continuous Adjustment Field of illumination : 70mm Filter : Cobalt Blue, Red Free, Heat Absorbing Filter Stand Section Type : Floor Stand Vertical Range Of Counterbalanced Arm : 400 mm General Information Base Size : 600mm X 550mm Weight : 50Kg. Height : 1550mm Electrical Power : AC 220V, 240V, 50/60 Hz Power Consumption : 220 VA/110V Optional Accessories Beam Splitter With CCD Camera Attachment CCD Camera Table Clamp Attachment Optional accessories are listed as separate items in "Accessories" section under products, there brought separately . Ent microscope, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,