



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fiber Optic Head Light

Product Code :
SAMCOMEDICAL4004A005



Description :

Fiber Optic Head Light Fiber Optic Head Light is adjustable and comfortable head light. It is incorporated with cushion that makes the the patient feel comfortable with the device. It is very durable and accurate in spot selection of the eye. Specifications Light Source Available with single/double outlet The front cable adaptor is suitable for CARL STORZ cable fitting For other makes adaptor is available separately (on request) Brilliant halogen illumination 24V-250W Compact and light weight instrument In built heat filter with extra white light with twin cooling fan inside Head Light Durable, comfortable head band fully adjustable Progressive IRISH diaphragm for spot selection Bright white light through front coated mirror without temperature increase Fiber Optic Cable 5mm dia CARL STORZ type adopter 7.5 feet length Fiber Optic Head Light, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,