



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Endo Camera (Digital Video Camera)

**Product Code :**  
SAMCOMEDICAL4004A006



**Description :**

Endo Camera (Digital Video Camera) It is very small and have a high source of light. It is highly flexible and has zooming features. It is lightweight and economic. Specifications Pick Up Device : 1/4" CCD Pixels : PAL = 522(H) X 582(V) Horizontal Resolution : NTSC/PAL = 450 Lines White Balance : Manuel White Balance Electrical Shutter : 1/60 (PAL = 50) - 1/96000 Power Source With 21 mm Indian Coupler : DC 12 V 500 (SMPS) Endo Camera (Digital Video Camera), ENT Equipment exporters, ENT Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),