



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Auto Projector Optical Equipment

**Product Code :**  
SAMCOMEDICAL4004A011



**Description :**

Auto Projector Optical Equipment ACP - 1800 Auto Projector Optical Equipment can provide 30 types of clear icons with different functions, and can program icon groups per specific needs. It can produce more precise measurement results in shorter time by remote controller and fast toggling icon, applying red/green & polarized icon test sheet, binocular eyesight and binocular balance. Features: Easy installation and maintenance It is quite easy to replace light bulbs, and unnecessary to adjust its focus and position. Customer Programmable Icons can be displayed in predefined sequence by programming. Display Single Icon Display single icon by multifunction mask plate. Tilted Placement Perfect image position can be obtained by adjusting the horizontal axis of projector. SPECIFICATION:- Projection distance 1.5m ~ 6m Projection magnification 30 x (at 5m) Projection size 330mm(W) x 270mm(H) (at 5m) Chart 30 different pattern Speed of chart conversion One chart per 0.03s Mask 1 open, 5 horizontal lines, 8 vertical lines, 21 single letters, 1 red/green speed of mask conversion One mask per 0.03s Program 2 sets programs, each program contains up to 30 steps Light source 12V, 50W Halogen Lamp Auto-off function After 10 minutes idle time. Power source AC 220V, 50Hz or 110V, 60Hz Power consumption 80W Dimensions 300mm(L) x 230mm(W) x 240mm(H) Net weight 6Kg Accessories Remote control, Polarized metal screen, halogen lamp, polarized glasses, fuses(2), batteries(2) Optional accessories Floor stand, wall stand Auto Projector Optical Equipment, ENT Equipment exporters, ENT Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

