



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Anti Bedsore Mattress

Product Code :
SAMCOMEDICAL10011A003

Description :

Anti Bedsore Mattress: Anti-Bedsore mattress (Air Jet Type). Tiny Holes on the surface for decreasing humidity. Replaces the foul air and reduce growth of microbes. It can effectively prevent from bedsore by dispersing weight of patient's body. Suitable for long term bedridden patients.



Samco Medical

Website: www.samcomedical.com,