



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Justice Auto Lensometer

Product Code :
SAMCOMEDICAL4004A013



Description :

Justice Auto Lensometer The only lensometer equipment at present that can measure accurately the low-degree lens below ± 0.12 in the country $\pm 1/2$ The lensometer equipment that can measure accurately the colored lens with a degree of transparency of 10% $\pm 1/2$ The lensometer equipment that can accurately and quickly measure the progressive lens. $\pm 1/2$ The only lensometer equipment that can display the pupil diameter $\pm 1/2$ Besides, the equipment can measure the pupil distance, and is equipped with UV meter and built-in thermosensitive printer $\pm 1/2$ It is characterized with fast measuring speed, clear interface, convenient operation, accuracy and good repeatability etc. Specifications: Sphere 0 to $\pm 25D$ (0.01D,0.12D,0.25D steps) Cylinder 0 to $\pm 10D$ (0.01D,0.12D,0.25D steps) Axis 0-180 $\pm 1/2$ (1 $\pm 1/2$ steps) ADD 0 to 10D (0.01D,0.12D,0.25D steps) Prism 0 to 10(0.01D,0.12, 0.25 steps) Step 0.01D PD 80mm(Optional) Lens Diameter Max.approx.120mm UV Out in UV (Optional) Compensation of Abbe number rang:20 to 60 Measurable >10% and over(20% and over $\pm 15D$ to $\pm 25D$) Contact lenses 0to $\pm 25D$ BC6.00 to 9.00 Display Full graphic LCD with back light(Black/White)320 ± 240 Printer Built-in thermosensitive printer, paper size 38mm(Optional) Power source 100-120V ~ /200-240V ~ choosing 50 or 60HZ Weight 7Kg Power AC220V/50HZ,AC110V/60Hz ,40W# Justice Auto Lensometer, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,