



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Auto Refractometer with Keratometer

Product Code :
SAMCOMEDICAL4004A015



Description :

Auto Refractometer with Keratometer We are engaged in offering Auto Refractometer with Keratometer. This equipment is used for quick and accurate keratometry measurement. Equipped SLD and highly sensitive CCD for high measurability of densely cataractous eyes and pseudophakic eyes, this equipment is designed using excellent quality components with latest technology. Features : $\frac{1}{2}$ With Auto-tracing feature, operator can adjust the focus and alignment in up and down direction $\frac{1}{2}$ With Auto video taping feature, operator can take measurements automatically once equipment obtains satisfactory focus and alignment $\frac{1}{2}$ With Servo Chin-Rest, operator can adjust height of Chin-rest easily $\frac{1}{2}$ Display is tilt-adjustable $\frac{1}{2}$ Printer has a paper auto-cut feature Specifications: $\frac{1}{2}$ Color screen/Auto test $\frac{1}{2}$ Vertex Distance(vd) 0.0, 12.0, 13.75mm $\frac{1}{2}$ Sphere Power(SPH) -20.00~+20.00D(in case of VD=12mm) $\frac{1}{2}$ Cylinder Power(CYL) 0.00~ $\frac{1}{2}$ 10.00D(0.12/0.25D Unit) $\frac{1}{2}$ Axis(AX) 1~1800 (10 Unit) $\frac{1}{2}$ Cylinder Form -, +, MIX $\frac{1}{2}$ Pupil Distance(PD) 10~80mm $\frac{1}{2}$ Minimum Pupil Diameter c2.00mm $\frac{1}{2}$ Radius of Curvature 5.0~10.2mm(Increments:0.01mm) Corneal Power 33.00~67.50D(When corner equivalent refractive Auto Refractometer with Keratometer, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,