



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Auto Refractometer Matronix Q 30 Koria

**Product Code :**  
SAMCOMEDICAL4004A016



**Description :**

Auto Refractometer Matronix Q 30 Koria We are engaged in offering Auto Refractometer Matronix q 30+ korean based technology. Technically advance with 3 year of warranty and money back garanty and has special IOL feature which enable to take the refractions of post operated patients. Features:  $\pm 1/2$  The most advanced technique of a graph edition, make the graph measured quickly and clear.  $\pm 1/2$  The precedent optical technique, consistent measurement and the accurate figure. Specifications: Refraction measurement Sphere (sph): -25D~ +22D, 0.12D/0.25D steps Cylinder (cyl): -80D~ +8D, 0.12D/0.25D steps Axis : 0~180 $\pm 1/2$ , 1 $\pm 1/2$  steps Pupil distance 45~85 mm, 0.5mm steps Vertical distance 0m, 12mm, 13.75mm Minimum pupil diameter 02.0mm Charter Automatic fogging method Display 5.7 inch true color LCD Print Import thermal printer Power Supply 100~120V/200~240V, 50/60HZ Dimension 275(W) \* 509(D)\*430-460(H)mm Weight 18Kg Auto Refractometer Matronix Q 30 Koria, ENT Equipment exporters, ENT Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),