



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
A BRAND NEW DIRECT OPHTHALMOSCOPES

Product Code :
SAMCOMEDICAL4004A027



Description :

A BRAND NEW DIRECT OPHTHALMOSCOPES OPHTHALMOSCOPES High quality optics for examining the interior structures of the eye, especially the retina, in a mini, pocket-size format OPHTHALMOSCOPES HEINE MINI 2000 LED critical tool offers all of the important features of a high-end ophthalmoscope combined in a compact, durable, portable configuration. Special Features: Enhanced XHL Xenon Halogen Technology; concentrated white light for perfect illumination & crisp retinal detail -20 to 20 with 18 lenses; -/ 1 to 1 1/4/ 4 in single steps 5 apertures for full diagnostic capability H 6.75 in. x W 1.25 in. x D 0.9 in. Weight: 3.7 oz A BRAND NEW DIRECT OPHTHALMOSCOPES, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,