



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
LENSOMETER LM- 510

**Product Code :**  
SAMCOMEDICAL4004A030



**Description :**

LENSOMETER LM- 510 Crossii nes with corona target External reading with tilling facility +/- to 25 Diopters power reading 0 to 5 prism reading 0 to 180" cylindrical axis Batteryandmainsoperated LENSOMETER LM- 510, ENT Equipment exporters, ENT Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),