



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
MOTORIZED TABLE (WITHOUT DRAWER)

**Product Code :**  
SAMCOMEDICAL4004A041



**Description :**

MOTORIZED TABLE (WITHOUT DRAWER) Dimension: 500(W)x500(D)x570(H) Weight: 29kg Highest: 880mm (i.e. 1/2 10mm) Lowest: 680mm (i.e. 1/2 10mm) Voltage: AC 110V or 220V (50/60Hz) Power: 75VA Fuse: 220V 5A (power) / 110V 7A (power) Table dimension: 480(W)x500(D)x24(T) MOTORIZED TABLE (WITHOUT DRAWER), ENT Equipment exporters, ENT Equipment manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),