



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
OT - TABLE

Product Code :
SAMCOMEDICAL4004A043



Description :

OT - TABLE This OT- Table being marketed by us is very advanced. It is a very heavy table made up of pure steel. The coverings are of very high quality materials. Besides, the table has sufficient space for the body to be kept in a comfortable position. This state of the art table can be fitted with all modern gadgets during the operation. And this is also not very expensive. Height Range min./max: 640 / 870 mm Input Definition: 230 V 50/60 Hz Lifting speed: 8 mm / sec 110 V 50/60 Hz Total length: 1880 mm Weight Capacity: 135 Kgs Total width: 690 mm Gross Weight: 115 Kgs Base size: 825 / 520 mm OT - TABLE, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,